

Brief Summary of the collaborative Project Entitled:
METABOLOMIC STUDY OF DISCORDANT PHENOTYPES OF
OBESITY AND TYPE 2 DIABETES

Biomedical Research Institute

Endocrinology and Nutrition Unit

Malaga Hospital Complex

Virgen de la Victoria Clinical Hospital

Málaga, Spain

Group Leader: Francisco Tinahones Madueño



Biomarkers and Nutritional & Food Metabolomics

Department of Nutrition and Food Science

Faculty of Pharmacy

University of Barcelona

Barcelona, Spain

Group Leader: Cristina Andrés-Lacueva


<http://www.nutrimetabolomics.com>

@nutrimetabolomics



The MHO Paradox

- ~ 25% Morbid Obese does NOT develop IR/T2D
- ~ 18% Not-Obese DOES develop IR/T2D



HPs - OBJs

The Blood and VAT Metabolome Contain Relevant Biological and Diagnostic Information

- To Identify Novel Biomarkers Able to Discriminate Discordant Phenotypes of Obesity/T2D from the 'Diabese' Phenotype and the Controls, to Unveil New Etiopathogenic mechanisms linking the two conditions

OBESITY

To be obese—does it matter if you are metabolically healthy?

Antony D. Karelis




nature
REVIEWS

ENDOCRINOLOGY

What Makes
The
Difference?

❖ WHAT DO WE HAVE ongoing? Study Designs & Samples Available (FIS-SAS)


Discordant Phenotypes: Observational STUDY

		pre T2D-T2D:	
		NO	YES
BMI (SEDO)		GROUP 1	GROUP 4 (DISC. PHENOTYPE A)
		GROUP 2	GROUP 5
		GROUP 3 (DISC. PHENOTYPE A)	GROUP 6

✓ Biosamples:


 n = 155

 n = 128

 n = 25 PILOT*

*3 groups; not exactly the same subjects

Gastric Bypass Follow-Up STUDY

		Time (months)			
		T0	T1	T3	T6
	pre T2D-	NO			
	T2D-	YES			

✓ n= 40 volunteers x 4 Timepoints
✓ Biosamples:

 n = 142
(some timepoint is missing)

Study Design

Observational STUDY

- ✓ Subjects recruited from Malaga Hospital Complex (Málaga, Spain); Informed consent; Ethical Committee approval
- ✓ INCLUSION Criteria / Grouping:

	Pre-diabetes or Diabetes: Hyperglycemia (fasting plasma glucose \geq mg/dL) Insulin Resistance (HOMA-IR $>$ 3.4)	
Body Mass Weight (BMI)	NO	YES
Normal Weight & Grade I Overweight (BMI $<$ 27)	GROUP 1	<u>GROUP 4</u> <u>(DISC. PHENOTYPE A)</u>
Grade I & II Obesity (BMI 30 - 40)	GROUP 2	GROUP 5
Morbid Obesity (BMI $>$ 40)	<u>GROUP 3</u> <u>(DISC. PHENOTYPE A)</u>	GROUP 6

- ✓ Sampling: n = aprox 30/group (sex- and age-matched)
- ✓ Biosamples: Blood Serum, Blood Plasma (pilot); Visceral Adipose Tissue
- ✓ EXCLUSION for: T2D medication, acute Inflammatory and Infective disease, alcohol & drug abuse

Potential Issues:

Treatment of confounders (MULTIdrugs, MULTIdiets)


Management of complexities (MULTIclass analysis by MVA)

Sampling: n $<$ in Discordant Phenotypes groups. Rare phenotypes



Discordant
Phenotypes:
Observational
STUDY

TARGETED APPROACH: Metabolic Profiling




≠ subjects
≠ samples

- ✓ QUANTITATIVE (concentrations nM – μM)
- ✓ **81** POLAR metabolites (Amino Acids, Biogenic Amines, Acylcarnitines, Total Hexoses)
- ✓ **399** NON-POLAR metabolites (Glycerophospholipids, Sphingolipids, Ceramides)
- ✓ Randomized sequence, QCs (30%)
- ✓ No replicates
- ✓ Metabolites filtered out if: a) CV>25% in QCs; b) >15% below LOD*

*McNewmar's test to compare the proportion of detectable values for discarded metabolites???

UNTARGETED APPROACH Metabolic Fingerprinting



- ✓ NOT QUANTITATIVE : mass signal intensities
- ✓ 5000 mass signals (m/z, RT) * ESI +/-
- ✓ Randomized sequence, QCs (30%)
- ✓ ISmix, ESmix. Apparent no need for normalization
- ✓ Low intensity signals:
- ✓ Noisy signals filtered out : excluded those that were NOT positively detected in at least the 40% of each of the study classes

High-throughput Transversal technologies -METABOLOMICS PLATFORM-

UNTARGETED APPROACH Metabolic Fingerprinting

- ✓ Hypothesis-GENERATING
- ✓ EXPLORATIVE, BIOMARKERS DISCOVERY
- ✓ MVA, Phenotype-Metabolome Associations
- ✓ IDEALLY unprocessed samples
- ✓ Exact Mass

LC-ESI-QToF-MS (QSTAR Elite)



TARGETED APPROACH: Metabolic Profiling

- ✓ Hypothesis-DRIVEN → only a cluster of metabolites
- ✓ QUANTITATIVE, BIOMARKERS VALIDATION
- ✓ UVA, ROC curves
- ✓ Pre-chromatographic treatment
- ✓ High sensitivity

LC-ESI-QqQ-MS/MS (MRM, PrIS, NL)



New QTRAP 6500
at CCIT-UB

<http://www.nutrimetabolomics.com>
candres@ub.edu